# KISS

Elektroden-Applikationssystem  
Electrode Application System  
Système d'application d'électrodes

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*A GE Medical Systems Company*

---

The product **Electrode Application System KISS** bears the CE mark

**CE-0366**

indicating its compliance with the provisions of the Council Directive 93/42/EEC concerning medical devices and fulfills the essential requirements of Annex I of this directive.

The product fully complies with the electromagnetic immunity requirements of standard IEC 601-1-2/EN 60601-1-2 "Electromagnetic Compatibility – Medical Electrical Equipment".

The device is in radio-interference protection class B in accordance with CISPR11/EN 55011.

The CE mark covers only the accessories listed in the "Equipment" section.

© Marquette Hellige GmbH 1998

Postfach 728

D-79007 Freiburg

Germany

Telephone +49 / 761 / 45 43-0

### Revision History

This document is subject to the Marquette Hellige change order system. The revision code, a letter that follows the document part number, changes with every update of the manual. The initial version of the manual is without revision code.

<b>Part No./Revision Code</b>	<b>Date</b>	<b>Comment</b>
227 450 01	27 August 19976	Initial Release
227 450 01-A	30 January 1998	ECO 203 198

## General Information

- \* This manual is an integral part of the device and describes its normal use. It should always be kept close to the device. Observance of the manual is a prerequisite for proper device performance and correct operation and ensures patient and operator safety.
- \* The symbol  denotes: Refer to Operator's Manual. It serves as an indicator for important facts to be noted when operating the device.
- \* Information which refers only to certain versions of the device is accompanied by the catalog number(s) of the instrument(s) concerned. The catalog number is given on the instrument nameplate.
- \* For detailed information on safety, please refer to the section "For your Safety" later on in this manual.
- \* Patient safety, specified measuring accuracy, and interference-free operation can be guaranteed only if original Marquette Hellige devices are interconnected (e.g. basic units and plug-in modules).
- \* Only use accessories which are listed in this user manual or which have been tested in combination with the device (e.g. patient cables, electrodes, transducers, sensors, consumables, etc.). If accessories or consumables from other manufacturers are used, Marquette Hellige cannot guarantee safe operation or functioning of the device.
- \* The warranty does not cover damages resulting from the use of accessories and consumables from other manufacturers.
- \* Marquette Hellige considers itself responsible for the effects on safety, reliability and performance of the device, only
  - if assembly operations, extensions, readjustments, modifications or repairs are carried out by Marquette Hellige or by persons authorized by Marquette Hellige
  - if the device is used in accordance with the instructions given in the user manual.
- \* For the address of the nearest Marquette Hellige sales-service office, please refer to the attached address list.
- \* All publications are in conformity with the device specifications and IEC publications on safety of electro-medical equipment valid at printing date. All rights are reserved for devices, circuits, techniques, software programs and names appearing in the manual.
- \* On request Marquette Hellige will supply a service manual.
- \* The Marquette Hellige management assurance system complies with the international DIN EN ISO 9001 and EN 46001 Standards.

# 1. Application and Functional Description

## 1.1 Introduction

The Marquette Hellige Electrode Application System KISS is a cord-integrated suction system and, as such, it is optimally suited to assist in the application of ECG electrodes for the recording of resting and stress-test ECGs. A pump which is either located in a separate housing or integrated in the EKG unit produces an adjustable partial vacuum in the electrodes, holding them in place on the skin at the respective electrode sites.

Via a common lead, the negative pressure and the electrical signals are conducted from the EKG unit to the electrode distributor. The electrode distributor may be equipped with 3, 10 or 12 electrode leads (see Table 1 below). It also accommodates the negative pressure control (60 to 220 mbar).

The electrodes are easily attached by pressing briefly on the suction cup. The suction pump starts up automatically when the first electrode is applied and continues to run until the adjusted negative pressure has been reached. After this, it operates only when necessary to correct a loss in pressure (caused, for instance, by the application of additional electrodes or by loose electrodes). When the negative pressure control is set at  $\odot$ , the pump switches off.

The ECG signal leads are equipped with overvoltage protection against defibrillation and HF surgery voltages.

The pump is powered from the power line (110 V to 240 V).

Electrode Connections	Electrodes	Leads
3 electrodes	R, L, N	EINTHOVEN, NEHB
10 electrodes	R, L, F, N, C1 to C6	EINTHOVEN + WILSON oder NEHB
12 electrodes	R, L, F, N, C1 to C6, Nax, Nst (C4=Nap)	EINTHOVEN + WILSON +NEHB

*Table 1*

The application system can be attached to an adjustable pole (accessory). Using special support systems, the pole can be mounted to the wall, to table tops or to equipment trolleys.


Figure 1-1. Electrode Application System KISS

## 1.2 For your Safety

- \* Dispose of the packaging material, observing the applicable waste control regulations and keeping it out of children's reach.
- \* The Application System **KISS** is designed to comply with IEC 601/EN 60601. It is a system of protection class II.
- \* Before connecting the device to the power line, check that the voltage and frequency ratings of the power line are similar to those indicated on the system's nameplate. If this is not the case, do not connect the system to the power line. Do not connect the system to extension cables with multiple power outlets.
- \* Before putting the system into operation visually inspect all connecting cables and electrodes for signs of damage. Damaged cables and electrodes must be replaced immediately.
- \* To disconnect the system from the power line, remove the plug from the wall outlet first.
- \* The ambient conditions defined in the Specifications section must be ensured.
- \* The Application System **KISS** is not suitable for operation in areas where an explosion hazard may occur.
- \* Magnetic and electrical fields are capable of interfering with the proper performance of the system. For this reason do not operate the Application System **KISS** in the vicinity of devices which generate strong, electromagnetic fields (X-ray or diathermia equipment or MRI devices).
- \* Do not connect the Application System **KISS** to devices that do not fulfill the requirements of standards IEC 601-1/EN 60601-1. Each of the single devices must be connected to a separate wall outlet. The use of electric power strips with multiple outlets is prohibited.
- \* Devices may only be interconnected with each other or with parts of systems when it has been made certain that there is no danger to the patient, the operators, or the environment as a result.  
In those instances where there is any element of doubt concerning the safety of connected devices, the user must contact the manufacturers concerned or other informed experts as to whether there is any possible danger to the patient, the operator, or the environment as a result of the proposed combination of devices. IEC standard 601-1 must be complied with in all cases.
- \* Liquids must not be allowed to enter the suction system or the pump housing. Components into which liquids have penetrated must be checked by a service technician before being used again.
- \* Use only the manufacturer's original accessories. Do not connect other signal sources to the cables.
- \* The classification into one of the categories B, BF or CF depends on the connected electrocardiograph classification. The same is true for the protection against defibrillation and HF surgery voltages (refer to the respective ECG recorder manuals). In general, extreme care must be exercised when applying medical electrical equipment. Many parts of the man/machine circuit are conductive, such as the patient, connectors, electrodes, transducers. It is very important that these conductive parts do not come into contact with other grounded, conductive parts when connected to the isolated patient input of the device. Such contact would bridge the patient's isolation and cancel the protection provided by the isolated input. In particular, there must be no contact of the neutral electrode and ground.
- \* Before using the system, the operator is obliged to ascertain that it is in correct working order and operating condition.
- \* The operator must be trained in the use of the equipment.
- \* Information which is of particular importance to the user is printed in a box throughout this manual.

## Literature

Medical Device Directive

IEC 601-1: 1990 + A1: 1993 + A2: 1995: Medical electrical equipment. General safety standards

IEC Publication 513/1994: Fundamental aspects of safety standards for medical electrical equipment.

## 2. Operating Controls and System Components


Figure 2-1. Application System KISS and pole

- |  |  |
|--|--|
| <ul style="list-style-type: none"> <li>1 Electrode distributor</li> <li>2 Negative pressure control</li> <li>3 Connector with electrode lead</li> <li>4 Silver/silver chloride electrode</li> <li>5 Filter disk</li> <li>6 Suction cup</li> <li>7 Electrocardiograph and pump connecting cable</li> <li>8 Connector for electrocardiograph</li> <li>9 Luer-Lock connector for suction pump</li> <li>10 Extension tube</li> <li>11 Power connector for pump module</li> <li>12 Pump module</li> <li>13 Screw clamp for fixation of pole to table top</li> <li>14 Wall-mount facility for pole</li> <li>15 Pole top with retaining screw, accommodates contact spray bottle</li> </ul> | <ul style="list-style-type: none"> <li>16 Locking catch to adjust the height of the swivel arm</li> <li>17 Locking catch to adjust the length of the swivel arm</li> <li>18 Bracket for electrode distributor</li> </ul> |
|--|--|

### Explanation of Symbols

-  Caution, refer to Operator's Manual
-  Protection class II
-  Function disabled
-  Adjustment of a variable parameter

### 3. System Setup

#### 3.1 Connecting the System


Figure 3-1. Connecting the application system to a *CardioSmart* with integrated suction pump

***The use of additional lengths of extension tubing is prohibited (overall tube length 4 meters/13 feet max.).***

When using an electrocardiograph or the *CORINA* model with integrated suction pump, all you have to do is plug in the patient cable connector (8) and connect the suction pump (9) (e.g. *CardioSmart*, Figure 3-1).

Should you be working with an application system whose pump is in a separate housing, proceed as follows:

- \* Find a suitable place for the pump module.
- \* Plug the patient cable connector (8) of the application system into the electrocardiograph.
- \* Connect the extension tube to connector (9) of the application system.
- \* Route the tube leading to the pump module, taking care that it is not kinked or squeezed.
- \* If necessary, cut the length of tube to size. Then slip the end of the tube onto the pump module nipple.
- \* Set the negative pressure control at the electrode distributor to the left stop (○) and connect the pump module to the power line. The module operates on line voltages from 110 to 240 V (switchover not necessary).

***With the negative pressure control at the left stop, the pump is switched off. The module itself remains in standby. In this mode the energy consumption is less than 0.5 W.***


Figure 3-2. Attaching the pole carrier to the wall-mount facility


Figure 3-3. Table mounting

## 3.2 Mounting the Pole

When mounting the pole, always take care to secure the pole top with the retaining screw either at the pole tube or at the pole carrier.

### Equipment Trolley

Please refer to the user manuals of the respective equipment trolleys for information on how to mount the pole to the various trolley models.

### Wall Mounting

If the application system is to be mounted on a wall, the pole carrier (a) must first be screwed on to the wall-mount facility (b). Then the wall-mount facility can be fixed at the desired location by means of four 8-mm screw inserts and the appropriate screws.

### Table Mounting

- \* Clamp the table mounting device to the table top as shown in Figure 3-3 or screw the top of the mounting device directly on to the table.
- \* The pole tube must be secured with a hexagon socket screw (arrow, Figure 3-3) (screw driver is included).


Figure 3-4. Putting the electrode distributor in place

### Mounting the Application System to the Pole

- \* Place the electrode distributor on the pins of the bracket (Figure 3-4).
- \* Route the joint suction and signal lead as shown in Figure 3-5.
- \* Secure the electrode leads as shown in Figure 3-5.
- \* When adjusting the height of the swivel arm, support the arm with one hand to prevent it from collapsing when you open clamp (16).


Figure 3-5. Securing the leads

**To avoid injuries which might be caused by the pins of the bracket, swing the swivel arm outside the work area when the electrode distributor is removed.**


Figure 4-1. Negative pressure control at electrode distributor  
a standard position

**☞ For standard applications the negative pressure control should be set to the gap in the arrow marker (a, Figure 4-1). This setting is in most cases sufficient for resting ECGs, while exercise tests will require a slightly higher setting. Be careful to determine the correct negative pressure level required in the respective situations. Be aware that excessive negative pressure settings and prolonged electrode application times may provoke hematomas (change electrode sites, if required). When using the system on small children or other patients with delicate skin, do not choose negative pressure settings in the upper range.**


Figure 4-2. Applying an electrode

**☞ In conjunction with the adapter 303 444 21, adhesive electrodes may be used instead of the suction electrodes.**

## 4. Operation

### 4.1 General Information

When connected to the power line, the pump module is automatically in standby mode. The slide control allows any negative pressure between 60 and 220 mbar to be adjusted. With the control set at  (left stop), the pump is switched off. The pump starts up automatically when a new electrode is applied, and it stops again as soon as the selected negative pressure has been reached. After this it will operate only to re-establish the partial vacuum required in the electrode. At the end of the ECG acquisition procedure, set the control to the  position. Air will then be automatically supplied to the system and the electrodes drop off. After about 3 seconds the pump switches off.

### 4.2 Applying the Electrodes

- \* Set the negative pressure control to the standard position (a, Figure 4-1).
- \* Apply a moderate amount of contact spray to the skin at the electrode site.
- \* Place the electrode on that site and press gently on the electrode's suction cup (Figure 4-2).

The pressure exerted on the cup causes the corresponding air valve to open, and the pump starts up. If the electrode has to be moved to another site, press on the sides of the suction cup to unseat the electrode: air will enter the cup and the air valve closes. Reapply the electrode to the new site by gently pressing on the suction cup. Always take care to select a negative pressure setting appropriate for the patient.

**☞ Use only small quantities of spray. Too much spray might prematurely clog the filter disk in the electrode. Do not spray the contact agent directly onto the electrodes. Do not use electrode cream or gel, as these substances may clog the air valve in the electrode distributor. Clean the electrodes as required and replace the filter disk as described below.**

 **At the end of the ECG acquisition procedure, always set the negative pressure control to the  position. This will extend the life of the pump.**

## 4.3 Ending the Procedure

- \* After the ECG acquisition procedure, set the negative pressure control to the left stop ( position).

As a result, air will enter the suction cups and the electrodes become detached. The pump switches off after 3 seconds.

Should you forget to set the control to the  position, the pump will switch off after 20 seconds. In this case you will have to set the control briefly to  before applying electrodes again.

- \* Clean the electrodes as described in section 6 “Cleaning and Maintenance”: **in case of heavy system usage, daily, otherwise every 3 to 7 days.**

### Replacing the Filter Disk

In case of heavy system usage, the filter disks should be replaced every day. Do not operate the system without filter disks!


Figure 4-3. Replacing the filter disk

- \* Disconnect the lead (c) from the electrode (d) by swiftly pulling off the lead connector, twisting it slightly.
- \* Remove the suction cup (a).
- \* Replace the filter disk (b) with a new one and reassemble the electrode.

## 5. System Test, Troubleshooting, Electrode Lead Replacement

### System Test

- \* Set the negative pressure control to a center position.
- \* Apply the electrode to be checked to your hand or arm.
- \* Set the negative pressure control to the right stop (max. setting): The dent of the suction cup must visibly enlarge.
- \* Move the negative pressure control further to the left to a weaker setting: The dent of the suction cup must visibly diminish.
- \* Set the negative pressure control to : Now it should be easy to lift off the electrode. The pump continues to operate another 3 seconds, then it stops.


You can conduct this test on each single electrode or on all electrodes in one pass.

For a test of the ECG signal transmission, the application system must be connected to an EKG unit.

## Troubleshooting

Problem	Problem	Remedy
Electrodes do not attach by suction or drop off, or the pump does not start up during electrode application	negative pressure control at 	set negative pressure control to the appropriate position
	pressure tube kinked or squeezed	route pressure tube correctly
	tube connecting the pump module interrupted	check connecting tube
	contaminated filter disk	replace filter disk
	interrupted power supply	power connector plugged in? supply voltage at wall outlet?
	contamination or moisture in air valve	replace electrode lead(s)
	faulty connection of electrode lead at distributor	check connections at distributor
	pump defect	notify sales-service office
The pump does not switch off when electrodes are detached	contamination or moisture in air valve	briefly set negative pressure control to  and start over
	electrode distributor not tight or contaminated	replace electrode lead notify sales-service office
12-Volt pump module: pump stops during procedure	power failure (depleted battery)	charge battery, check connections

*Figure 5-1. Replacing an electrode lead*

### Replacing Electrode Leads

The electrode leads are subject to wear and tear and, therefore, they can be replaced. They are supplied including the air valve, but excluding the electrode (section 8 “Equipment”). The air valve is integrated in the connector that is plugged into the distributor.

- \* Remove the electrode from the damaged lead with a slight twist.
- \* Remove the damaged lead from the electrode distributor, grasping the lead directly at the connection (Figure 5-1).
- \* Plug the connector of the new lead into the electrode distributor.

Always hold the connector as shown in Figure 5-1. Do not force the connector into the socket and take care that it is properly inserted. Otherwise there may be problems with the tightness of the system and with the necessary contacts. Also take care to insert the leads in the correct sockets, observing the lead designations printed on the distributor housing.

- \* Connect the electrode to the new lead.

## 6. Cleaning and Maintenance

### 6.1 Cleaning and Disinfection

#### Instrument Surface

- \* Disconnect the system from the power line before cleaning and disinfecting the instrument surface (pump module).
- \* Liquids must not be allowed to enter the module. Use a moist cloth to wipe the module clean. You may use any hospital-grade cleaning agents and disinfectants.

**Disinfectants with a phenol or peroxide base are unsuitable for disinfection of the module surface.**

#### Leads

- \* Disconnect the cables from the system components before cleaning and disinfecting them. To do so, always pull on the plug and not on the cable.
- \* Rub the cables clean with a cloth soaked in soapy water. Disinfect the cables by rubbing them down with disinfectant.

**Under no circumstances should cables be immersed in the cleaning solution or disinfectant. Even if the system were not under negative pressure, the liquid would enter the system, spoiling the filter disks and clogging the electrode application system.**

#### Electrodes

- \* In case of heavy usage of the system, clean the electrodes every day, otherwise clean them every 3 to 7 days. Do not forget to replace the filter disks at the same time (part no. 932 047 75, 100/cs).
- \* The surface of electrodes with minor contamination can be cleaned with a cloth soaked in water and detergent. To disinfect the electrodes, wipe them down with a cloth moistened in 80% alcohol.

**Do not immerse the electrodes in liquids.**

- \* For thorough cleaning, disconnect the electrode (d) from the lead (c) with a slight twist.
- \* Next, remove the suction cup (a) and the filter disk (b) (Figure 6-1).


Figure 6-1. Cleaning an electrode

- \* Clean the electrode with water and a detergent. Use a small brush (toothbrush) to remove grime.

**Do not use pointed metal objects to remove solid particles of dirt from the electrodes as this would destroy the silver/silver-chloride layer. Use a commercially available fiberglass eraser instead.**

- \* All commercially available cleaning agents used for surgical instruments are suitable for cleaning of the electrodes. Follow the manufacturer's instructions to mix the preparation. Do not use metal dishes and take care not to immerse plugs and metal sockets in the solution.
- \* Carefully dry the electrodes. Then reassemble them, fitting a new filter disk.

### 6.2 Maintenance

The pump module does not require any maintenance.

## 7. Specifications

Electrode application system for acquiring ECG signals from a patient

### Electrode distributor

for 3, 10 or 12 plug-in electrode leads

### Electrodes

silver-silver chloride electrodes, with connectors

### Defibrillation protection

by means of integrated current and voltage limitation

### Signal connector

15-pin connection, pin configuration

- 1 = C2
- 2 = C3
- 3 = C4 = Nap
- 4 = C5
- 5 = C6
- 6 = shield
- 7 = bridge to 8 (only with 12 electrodes)
- 8 = bridge to 7 (only with 12 electrodes)
- 9 = R
- 10 = L
- 11 = F
- 12 = C1
- 13 = Nst
- 14 = N
- 15 = Nax

### Negative pressure adjustment for pump module

continuously adjustable between 60 and 220 mbar;  
in position  negative pressure < 30 mbar

### Power supply

from the power line; design in compliance with IEC 601/EN 60601, protection class II;

for rated line voltages from 110 to 240 V, 49 to 65 Hz (without switchover)

operating voltage range 98 to 264 V, AC  
rated current range 5 to 100 mA,  
in standby mode approx. 20 mA

### Mechanical data

swivel arm length (extendible) 75 to 125 cm/29.5 to 49 in.  
vertical swivel range 45°  
horizontal swivel range 360°

length of patient cable from  
electrode distributor  
to EKG unit approx. 250 cm/8.2 ft.

length of electrode leads limbs 140 cm/4.6 ft.  
chest 120 cm/3.9 ft.

length of negative pressure tube  
connected to separate pump 200 cm/6.6 ft.

### Environment

#### *Operation*

- \* temperature between +10 to +40 °C/+50 to +104 °F
- \* relative humidity between 25 and 95%, non-condensing
- \* atmospheric pressure between 700 and 1060 hPa

#### *Storage and transport*

- \* temperature between -30 to +60 °C/-22 to +140 °F
- \* relative humidity between 25 and 95%, non-condensing
- \* atmospheric pressure between 500 and 1060 hPa

**Dimensions** (pump module)

- \* length 130 mm/5.1 in.
- \* width 65 mm/2.6 in.
- \* height 70 mm/2.8 in.

**Weight**

- \* electrode distributor  
with 12 electrodes approx. 800 g/1.7 lb.
- pump module 700 g/1.5 lb.

## 8. Equipment

Subject to modification.

Always refer to latest list of accessories.

216 120 ..	Electrode Application System KISS, 3 electrode connections, w/o. pump	303 443 65	Electrode lead C1
		303 443 66	Electrode lead C2
216 121 ..	Electrode Application System KISS, 10 electrode connections, w/o. pump	303 443 67	Electrode lead C3
		303 443 68	Electrode lead C4
216 122 ..	Electrode Application System KISS, 12 electrode connections, w/o. pump	303 443 69	Electrode lead C5
		303 443 70	Electrode lead C6
216 123 ..	Electrode Application System KISS, 3 electrode connections, with pump	303 443 71	Electrode lead F
		303 443 72	Electrode lead L
216 124 ..	Electrode Application System KISS, 10 electrode connections, with pump	303 443 73	Electrode lead R
		303 443 74	Electrode lead N
216 125 ..	Electrode Application System KISS, 12 electrode connections, with pump	303 443 75	Electrode lead N <sub>st</sub>
227 450 01	Operator's Manual	303 443 76	Electrode lead N <sub>ax</sub>
303 443 77	Swivel arm	303 444 21	Clip adapter for adhesive electrodes
384 015 84	Table mount system with pole tube	303 444 19	Extension tube for connection of pump
384 013 30	Wall mount system for swivel arm		
303 444 20	Ergometer mount system		
217 307 01	Electrode contact spray, 200-ml bottle		
217 307 05	Electrode contact spray, refill bottle		
932 047 75	Filter disks (100(cs.))		
217 328 01	Suction electrode, assembled		
432 521 25	Suction cup		

## EC Declaration of Conformity

Document No. 33-95

**HELLIGE GMBH, Munzinger Strasse 3, D-79111 Freiburg, Germany**

We herewith declare that the product

**ECG Electrode Application System ES 600 (KISS x), product status Version 1.0**

fulfills the requirements of the following directives, standards and normative documents:

1. Council Directive 93/42/EEC of 14 June 1993 concerning medical devices
2. EN 60601-1:1990 / DIN VDE 0750 Teil 1/12.91 / IEC 601-1 2. Ausgabe 1988  
EN 60601-1-2:1993 / DIN EN 60601-1-2:1994-09 / IEC 601-1-2:1993  
IEC 601-2-25:1993 / E DIN IEC 601-2-25 (Klassifikation VDE 0750 Teil 2-25):1994-11

Compliance of a representative sample of the designated product with the "essential requirements" of Annex I of the Directive 93/42/EEC has been certified by

**HELLIGE GMBH, Quality Management and Certification,  
Munzinger Strasse 3, D-79111 Freiburg, Germany  
Test Report No. PB95002 of 21 December 1995.**

The medical device has been assigned to class **IIb** as specified in Annex IX of the Directive 93/42/EEC. It bears the marking

**CE-0366.**

The designated product has been designed, manufactured and tested under a quality management system according to DIN ISO 9001 / EN 29001 and EN 46001.  
The conformity of the quality management system has been certified by:

**VDE Testing and Certification Institute**

Freiburg, 21 December 1995

HELLIGE GMBH

Karl Braun, Director  
Marketing & Engineering

The technical documentation  
is filed in the QA Dept.

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